

## VITA VM9 Classical Powder Indication Table

VITAPAN Classical Shades A1-D4	Coloring Liquid	Effect Bonder	Margin	Effect Liner	Chroma Plus	Enamel
A1	LL1	EB1	M1/M7*	EL2	CP1	ENL
A2	LL2	EB2	M4/M7*	EL1/EL3*	CP2	ENL
A3	LL2	EB2	M4	EL4/EL6*	CP2/CP3*	ENL
A3.5	LL3	EB3	M4/M9*	EL5/EL6*	CP2/CP3*	END
A4	LL3	EB3	M4/M9*	EL1/EL4*	CP2/CP4*	END
B1	LL1	EB1	M1/M4*	EL1/EL2*	CP1	END
B2	LL1	EB1	M1/M4*	EL1/EL2*	CP1	END
B3	LL3	EB3	M4	EL2/EL4*	CP2/CP3*	END
B4	LL3	EB3	M4/M9*	EL4/EL6*	CP3	END
C1	LL3	EB3	M1/M4*	EL1/EL6*	CP1	END
C2	LL2	EB2	M4/M7*	EL2/EL6*	CP1/CP5*	END
C3	LL3	EB3	M4/M7*	EL6	CP1/CP5*	ENL
C4	LL4	EB4	M4/M7*	EL3/EL6*	CP5	ENL
D2	LL2	EB2	M1/M9*	EL2/EL6*	CP1/CP5*	END
D3	LL3	EB3	M4/M7*	EL2/EL3*	CP2/CP5*	END
D4	LL3	EB3	M1/M4*	EL2/EL6*	CP2/CP5*	END

\* Mixing Ratio 1:1

