

# VITAVM<sup>®</sup>13 porcelain now available in VITA 3D-Master<sup>®</sup> and VITA Classical shades



## VITAVM<sup>®</sup>13

**VITA<sup>®</sup>**

VITA's popular VM13 veneering porcelain offers outstanding physical properties, handling characteristics and, due to its fine microstructure, clinical wear characteristics that mimic those of enamel.

- Low firing temperature and reduced shrinkage
- Excellent shade coverage with high chroma, high opacity Base Dentines
- CTE of 13.1-13.6 for veneering metal frameworks with a CTE of 13.8-15.2
- Extremely low antagonist tooth wear rates
- Optimum firing stability, which preserves shape and shade even after several firings

**800-828-3839**

800-263-4778 in Canada • [www.vident.com](http://www.vident.com)

**VIDENT<sup>™</sup>**

Classical Kits

### VM 13 Basic Kit Classical Powder V50-0001

- 1 VITA VM 13 Wash Opaque
- 16 VITA VM 13 Opaque
- 16 VITA VM 13 Base Dentine
- 5 VITA VM 13 Chroma Plus
- 2 VITA VM 13 Enamel
- 1 VITA VM 13 Neutral
- 1 VITA VM 13 Window
- 3 VITA VM 13 Corrective
- 1 VITA VM Modeling Liquid, 50 ml
- 1 VITA VM Opaque Fluid, 50 ml
- Auxiliary Products

### VM 13 Basic Kit Classical Paste V50-0002

- 1 VITA VM 13 Wash Opaque Paste
- 16 VITA VM 13 Opaque Paste
- 16 VITA VM 13 Base Dentine
- 5 Chroma Plus
- 2 VITA VM 13 Enamel
- 1 VITA VM 13 Neutral
- 1 VITA VM 13 WIN
- 3 VITA VM 13 Corrective
- 1 VITA VM Modeling Liquid, 50 ml
- 1 VITA VM Paste Fluid, 15 ml
- Auxiliary Products

### VM 13 Basic Kit Classical Powder (Small) V50-0003

- 1 VITA VM 13 Wash Opaque
- 6 VITA VM 13 Opaque\*
- 6 VITA VM 13 Base Dentine\*
- 5 VITA VM 13 Chroma Plus
- 2 VITA VM 13 Enamel
- 1 VITA VM 13 Neutral
- 1 VITA VM 13 Window
- 1 VITA VM 13 Corrective
- 1 VITA VM Modeling Liquid, 50 ml
- 1 VITA VM Opaque Fluid, 50 ml
- Auxiliary Products

### VM 13 Basic Kit Classical Paste (Small) V50-0004

- 1 VITA VM 13 Wash Opaque Paste
- 6 VITA VM 13 Opaque Paste\*
- 6 VITA VM 13 Base Dentine\*
- 5 Chroma Plus
- 2 VITA VM 13 Enamel
- 1 VITA VM 13 Neutral
- 1 VITA VM 13 Window
- 1 VITA VM 13 Corrective
- 1 VITA VM Modeling Liquid, 50 ml
- 1 VITA VM Paste Fluid, 15 ml
- Auxiliary Products

### VM 13 Classical Buildup Kit V50-0005

- 16 VITA VM 13 Transpa Dentine
- 1 VITA VM Modeling Liquid, 50 ml

### VM 13 Classical Buildup Kit (Small) V50-0006

- 6 VITA VM 13 Transpa Dentine
- 1 VITA VM Modeling Liquid, 50 ml

### VM 13 Classical Color Kit Paste V50-0007

- 16 VITA VM 13 Opaque Paste
- 16 VITA VM 13 Base Dentine
- 16 VITA VM 13 Transpa Dentine
- 1 VITA VM Modeling Liquid, 50ml
- 1 VITA VM Paste Fluid, 15 ml
- 1 Classical Shade Guide
- 1 VITA VM13 Shade Indicator Classical

### VM 13 Classical Color Kit, Powder V50-0008

- 16 VITA VM13 Classical Opaque
- 16 VITA VM13 Classical Base Dentine
- 16 VITA VM13 Classical Transpa Dentine
- 1 VITA VM Modeling Liquid, 50ml
- 1 VITA VM Opaque Fluid, 50ml
- 1 Classical Shade Guide
- 1 VITA VM13 Shade Indicator

### VM 13 One Color Kit, A2 Powder V50-0011

- 1 VITA VM 13 Wash Opaque
- 1 VITA VM13 Opaque
- 1 VITA VM 13 Base Dentine
- 1 VITA VM 13 Transpa Dentine
- 1 VITA VM 13 Enamel
- 1 VITA VM 13 Neutral
- 1 VITA VM 13 Window
- 1 VITA VM Modeling Liquid, 50ml
- 1 VITA VM Opaque Fluid, 50ml

### VM 13 One Color Kit, A3 Powder V50-0012

- 1 VITA VM 13 Wash Opaque
- 1 VITA VM13 Opaque
- 1 VITA VM 13 Base Dentine
- 1 VITA VM 13 Transpa Dentine
- 1 VITA VM 13 Enamel
- 1 VITA VM 13 Neutral
- 1 VITA VM 13 Window
- 1 VITA VM Modeling Liquid, 50ml
- 1 VITA VM Opaque Fluid, 50ml

\* A1, A2, A3, A3.5, B3, D3

3D-Master Kits

### VM 13 Basic Kit Powder Opaque VX50-0001

- 1 VITA VM 13 Wash Opaque Powder
- 5 VITA VM 13 Opaque Powder
- 26 VITA VM 13 Base Dentine Powder
- 5 VITA VM 13 Chroma Plus Powder
- 2 VITA VM 13 Enamel Powder
- 1 VITA VM 13 Neutral Powder
- 1 VITA VM 13 Window Powder WIN
- 3 VITA VM 13 Corrective Powder
- 1 VITA VM 13 Opaque Fluid, 50ml
- 1 Modeling Liquid, 50ml
- Auxiliary Products

### VM 13 Basic Kit Paste Opaque VX50-0002

- 1 VITA VM 13 Wash Opaque Paste
- 5 VITA VM 13 Paste Opaque
- 26 VITA VM 13 Base Dentine Powder
- 5 VITA VM 13 Chroma Plus Powder
- 2 VITA VM 13 Enamel Powder
- 1 VITA VM 13 Neutral Powder
- 1 VITA VM 13 Window Powder WIN

- 3 VITA VM 13 Corrective Powder
- 1 VITA VM 13 Paste Fluid, 15ml
- 1 Modeling Liquid, 50ml
- Auxiliary Products

### VM 13 Basic Kit Powder Opaque (Small) VX50-0003

- 1 VITA VM 13 Wash Opaque Powder
- 4 VITA VM 13 Opaque Powder
- 15 VITA VM 13 Base Dentine Powder
- 5 VITA VM 13 Chroma Plus Powder
- 2 VITA VM 13 Enamel Powder
- 1 VITA VM 13 Neutral Powder
- 1 VITA VM 13 Window Powder WIN
- 1 VITA VM 13 Corrective Powder
- 1 VITA VM 13 Opaque Fluid, 50ml
- 1 Modeling Liquid, 50ml
- Auxiliary Products

### VM 13 Basic Kit Paste Opaque (Small) VX50-0004

- 1 VITA VM 13 Wash Opaque Paste
- 4 VITA VM 13 Paste Opaque

- 15 VITA VM 13 Base Dentine Powder
- 5 VITA VM 13 Chroma Plus Powder
- 2 VITA VM 13 Enamel Powder
- 1 VITA VM 13 Neutral Powder
- 1 VITA VM 13 Window Powder
- 1 VITA VM 13 Corrective Powder
- 1 VITA VM 13 Paste Fluid, 15ml
- 1 Modeling Liquid, 50ml
- Auxiliary Products

### VM 13 Buildup Kit VX50-0005

- 26 VITA VM 13 Transpa Dentine
- 1 Modeling Liquid, 50ml
- 1 Working Instructions

### VM 13 Buildup Kit (Small) VX50-0006

- 15 VITA VM 13 Transpa Dentine
- 1 Modeling Liquid, 50ml
- 1 Working Instructions

### VM 13 Bleached Shade Kit Powder Opaque VX50-0007

- 1 VITA VM 13 Opaque Powder
- 3 VITA VM 13 Base Dentine Powder
- 3 VITA VM 13 Transpa Dentine Powder
- 1 VITA VM 13 Enamel Light ENL
- 1 VITA VM 13 Neutral Powder NT
- 1 VITA VM 13 Window Powder
- 1 Modeling Liquid, 50ml
- 1 Opaque Fluid, 50ml
- Auxiliary Products

### VM 13 Bleached Shade Kit Paste Opaque VX50-0008

- 1 VITA VM 13 Paste Opaque
- 3 VITA VM 13 Base Dentine Powder
- 3 VITA VM 13 Transpa Dentine Powder
- 1 VITA VM 13 Enamel Light
- 1 VITA VM 13 Neutral Powder
- 1 VITA VM 13 Window Powder
- 1 Modeling Liquid, 50ml
- 1 VITA VM 13 Paste Fluid, 15ml
- Auxiliary Products

Specialty Kits

### VM 13 Professional Kit VX50-0009

- 11 VITA VM 13 Effect Chroma
- 11 VITA VM 13 Effect Enamel
- 3 VITA VM 13 Mamelon
- 3 VITA VM 13 Effect Pearl
- 3 VITA VM 13 Effect Opal
- 6 VITA VM 13 Effect Liner
- 3 VITA VM13 Color Indicators

### VM 13 Professional Kit (Small) VX50-0010

- 5 VITA VM 13 Effect Chroma
- 1 VITA VM 13 Mamelon
- 1 VITA VM 13 Effect Pearl
- 1 VITA VM 13 Effect Opal
- 7 VITA VM 13 Effect Enamel
- 1 VITA VM13 Color Indicator

### VM 13 Gingiva Kit Powder Opaque VX50-0017

- 5 VITA VM 13 Gingiva
- 2 VITA VM 13 Gingiva Opaque
- 1 VITA VM13 Color Indicator

### VM 13 Gingiva Kit Paste Opaque VX50-0018

- 5 VITA VM 13 Gingiva
- 2 VITA VM 13 Gingiva Paste Opaque
- 1 VITA VM 13 Color Indicator

### VM 13 Margin Kit VX50-0019

- 6 VITA VM 13 Margin
- 1 Color Indicator

**800-828-3839 • 800-263-4778 in Canada • www.vident.com**